

Product Datasheet

Bovine TTR ELISA Kit (orb1197040)

Description	Bovine Transthyretin (TTR) ELISA Kit
Reactivity	Bovine
Range	7.81-500ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	competitive inhibition
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O46375
Sensitivity	2.8ng/mL
Entrez	280948
Expiration Date	Please enquire.