

Product Datasheet

Bovine Tie2 ELISA Kit (orb1197057)

Description	Bovine TEK Tyrosine Kinase, Endothelial (Tie2) ELISA Kit
Reactivity	Bovine
Range	6.25-400ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q06807
Sensitivity	1.56ng/mL
Entrez	280939
Expiration Date	Please enquire.