

Product Datasheet

Bovine PPARg ELISA Kit (orb1197088)

Description	Bovine Peroxisome Proliferator Activated Receptor Gamma (PPARg) ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O18971
Sensitivity	0.115ng/mL
Entrez	281993
Expiration Date	Please enquire.