

Product Datasheet

Bovine GAL9 ELISA Kit (orb1197164)

Description	Bovine Galectin 9 (GAL9) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q3MHZ8
Sensitivity	7.8pg/mL
Entrez	510813
Expiration Date	Please enquire.