

Product Datasheet

Bovine GAL8 ELISA Kit (orb1197165)

Description	Bovine Galectin 8 (GAL8) ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q2YDD7
Sensitivity	0.074ng/mL
Entrez	506977
Expiration Date	Please enquire.