

Product Datasheet

Bovine GAL4 ELISA Kit (orb1197167)

Description	Bovine Galectin 4 (GAL4) ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q3T0D6
Sensitivity	0.112ng/mL
Entrez	614804
Expiration Date	Please enquire.