

Product Datasheet

Porcine FAS ELISA Kit (orb1197172)

Description	Porcine Factor Related Apoptosis (FAS) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	077736
Sensitivity	7.2pg/mL
Entrez	396826
Expiration Date	Please enquire.