

Product Datasheet

Bovine DEFb2 ELISA Kit (orb1197176)

Description	Bovine Defensin Beta 2 (DEFb2) ELISA Kit
Reactivity	Bovine
Range	62.5-4000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P46160
Sensitivity	25.7pg/mL
Entrez	404063
Expiration Date	Please enquire.