

Product Datasheet

Bovine CYPB ELISA Kit (orb1197178)

Description	Bovine Cyclophilin B (CYPB) ELISA Kit
Reactivity	Bovine
Range	12.5-800ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P80311
Sensitivity	5.2ng/mL
Entrez	281419
Expiration Date	Please enquire.