

Product Datasheet

Porcine CXCR4 ELISA Kit (orb1197179)

Description	Porcine Chemokine C-X-C-Motif Receptor 4 (CXCR4) ELISA Kit
Reactivity	Porcine
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q764M9
Sensitivity	0.059ng/mL
Entrez	396659
Expiration Date	Please enquire.