

Product Datasheet

Porcine CX43 ELISA Kit (orb1197180)

Description	Porcine Connexin 43 (CX43) ELISA Kit
Reactivity	Porcine
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9GJX4
Sensitivity	0.16ng/mL
Entrez	100518636
Expiration Date	Please enquire.