

Product Datasheet

Bovine CTSL ELISA Kit (orb1197182)

Description	Bovine Cathepsin L (CTSL) ELISA Kit
Reactivity	Bovine
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	A4IFS7
Sensitivity	0.314ng/mL
Entrez	515200
Expiration Date	Please enquire.