

Product Datasheet

Bovine CRP ELISA Kit (orb1197186)

Description	Bovine C Reactive Protein (CRP) ELISA Kit
Reactivity	Bovine
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	C4T8B4
Sensitivity	6.2pg/mL
Entrez	527553
Expiration Date	Please enquire.