

Product Datasheet

Bovine CKM ELISA Kit (orb1197190)

Description	Bovine Creatine Kinase, Muscle (CKM) ELISA Kit
Reactivity	Bovine
Range	12.5-800U/L
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9XSC6
Sensitivity	4.65U/L
Entrez	286822
Expiration Date	Please enquire.