

Product Datasheet

Porcine CDHE ELISA Kit (orb1197191)

Description	Porcine Cadherin, Epithelial (CDHE) ELISA Kit
Reactivity	Porcine
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	C6EVT4
Sensitivity	0.195ng/mL
Entrez	100048953
Expiration Date	Please enquire.