

Product Datasheet

Bovine CDHE ELISA Kit (orb1197192)

Description	Bovine Cadherin, Epithelial (CDHE) ELISA Kit
Reactivity	Bovine
Range	0.937-60ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q6R8F2
Sensitivity	0.23ng/mL
Entrez	282637
Expiration Date	Please enquire.