

Product Datasheet

Bovine APOH ELISA Kit (orb1197209)

Description	Bovine Apolipoprotein H (APOH) ELISA Kit
Reactivity	Bovine
Range	3.12-200ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P17690
Sensitivity	1.87ng/mL
Entrez	281006
Expiration Date	Please enquire.