

Product Datasheet

Bovine ANGPT1 ELISA Kit (orb1197219)

Description	Bovine Angiopoietin 1 (ANGPT1) ELISA Kit
Reactivity	Bovine
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	O18920
Sensitivity	0.195ng/mL
Entrez	282140
Expiration Date	Please enquire.