

Product Datasheet

Bovine AMH ELISA Kit (orb1197221)

Description	Bovine Anti-Mullerian Hormone (AMH) ELISA Kit
Reactivity	Bovine
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P03972
Sensitivity	0.117ng/mL
Entrez	280718
Expiration Date	Please enquire.