

Product Datasheet

Porcine FGF21 ELISA Kit (orb1197240)

Description	Porcine Fibroblast Growth Factor 21 (FGF21) ELISA Kit
Reactivity	Porcine
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	C6JWG4
Sensitivity	13.3pg/mL
Entrez	100302504
Expiration Date	Please enquire.