

## Product Datasheet

### Rat SPJ ELISA Kit (orb1197259)

|                            |                                                     |
|----------------------------|-----------------------------------------------------|
| <b>Description</b>         | Rat Surfactant Associated Protein J (SPJ) ELISA Kit |
| <b>Reactivity</b>          | Rat                                                 |
| <b>Range</b>               | 0.156-10ng/mL                                       |
| <b>Tested Applications</b> | ELISA                                               |
| <b>Note</b>                | For research use only                               |
| <b>Application notes</b>   | sandwich                                            |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids.      |
| <b>Uniprot ID</b>          | <b>Q5I0E5</b>                                       |
| <b>Sensitivity</b>         | 0.061ng/mL                                          |
| <b>Entrez</b>              | <b>315953</b>                                       |
| <b>Expiration Date</b>     | Please enquire.                                     |