

Product Datasheet

Rat SEPP1 ELISA Kit (orb1197268)

Description	Rat Selenoprotein P1, Plasma (SEPP1) ELISA Kit
Reactivity	Rat
Range	6.25-400ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P25236
Sensitivity	2.24ng/mL
Entrez	29360
Expiration Date	Please enquire.