

Product Datasheet

Rat PRF1 ELISA Kit (orb1197281)

Description	Rat Perforin 1 (PRF1) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q5FVS5
Sensitivity	0.18ng/mL
Entrez	50669
Expiration Date	Please enquire.