

Product Datasheet

Rat MSTN ELISA Kit (orb1197312)

Description	Rat Myostatin (MSTN) ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	035312
Sensitivity	0.214ng/mL
Entrez	29152
Expiration Date	Please enquire.