

Product Datasheet

Rat MGP ELISA Kit (orb1197317)

Description	Rat Matrix Gla Protein (MGP) ELISA Kit
Reactivity	Rat
Range	93.7-6000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P08494
Sensitivity	39.65pg/mL
Entrez	25333
Expiration Date	Please enquire.