

Product Datasheet

Mouse GC ELISA Kit (orb1197358)

Description	Mouse Glucagon (GC) ELISA Kit
Reactivity	Mouse
Range	23.4-1500pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P55095
Sensitivity	17.34pg/mL
Entrez	14526
Expiration Date	Please enquire.