

Product Datasheet

Rat COMP ELISA Kit (orb1197393)

Description	Rat Cartilage Oligomeric Matrix Protein (COMP) ELISA Kit
Reactivity	Rat
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P35444
Sensitivity	0.65ng/mL
Entrez	25304
Expiration Date	Please enquire.