

Product Datasheet

Rat CBG ELISA Kit (orb1197405)

Description	Rat Corticosteroid Binding Globulin (CBG) ELISA Kit
Reactivity	Rat
Range	3.12-200ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P31211
Sensitivity	1.39ng/mL
Entrez	299270
Expiration Date	Please enquire.