

## Product Datasheet

### Rat CACT ELISA Kit (orb1197410)

|                            |                                                          |
|----------------------------|----------------------------------------------------------|
| <b>Description</b>         | Rat Carnitine Acylcarnitine Translocase (CACT) ELISA Kit |
| <b>Reactivity</b>          | Rat                                                      |
| <b>Range</b>               | 0.312-20ng/mL                                            |
| <b>Tested Applications</b> | ELISA                                                    |
| <b>Note</b>                | For research use only                                    |
| <b>Application notes</b>   | sandwich                                                 |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids.           |
| <b>Uniprot ID</b>          | <b>P97521</b>                                            |
| <b>Sensitivity</b>         | 0.087ng/mL                                               |
| <b>Entrez</b>              | <b>117035</b>                                            |
| <b>Expiration Date</b>     | Please enquire.                                          |