

## Product Datasheet

### Rat C5 ELISA Kit (orb1197411)

|                            |                                              |
|----------------------------|----------------------------------------------|
| <b>Description</b>         | Rat Complement Component 5 (C5)<br>ELISA Kit |
| <b>Reactivity</b>          | Rat                                          |
| <b>Range</b>               | 0.312-20ng/mL                                |
| <b>Tested Applications</b> | ELISA                                        |
| <b>Note</b>                | For research use only                        |
| <b>Application notes</b>   | sandwich                                     |
| <b>Sample Types</b>        | serum, plasma or other biological fluids.    |
| <b>Uniprot ID</b>          | <b>P08650</b>                                |
| <b>Sensitivity</b>         | 0.138ng/mL                                   |
| <b>Entrez</b>              | <b>362119</b>                                |
| <b>Expiration Date</b>     | Please enquire.                              |