

Product Datasheet

Rat BRAK ELISA Kit (orb1197412)

Description	Rat Breast And Kidney Expressed Chemokine (BRAK) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q8K453
Sensitivity	0.058ng/mL
Entrez	306748
Expiration Date	Please enquire.