

Product Datasheet

Rat AST ELISA Kit (orb1197416)

Description	Rat Aspartate Aminotransferase (AST) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P13221
Sensitivity	0.117ng/mL
Entrez	24401
Expiration Date	Please enquire.