

Product Datasheet

Rat APAF1 ELISA Kit (orb1197425)

Description	Rat Apoptotic Peptidase Activating Factor 1 (APAF1) ELISA Kit
Reactivity	Rat
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9EPV5
Sensitivity	0.062ng/mL
Entrez	78963
Expiration Date	Please enquire.