

Product Datasheet

Rat ALOX5 ELISA Kit (orb1197435)

Description	Rat Arachidonate-5-Lipoxygenase (ALOX5) ELISA Kit
Reactivity	Rat
Range	0.781-50ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P12527
Sensitivity	0.35ng/mL
Entrez	25290
Expiration Date	Please enquire.