

## Product Datasheet

### Rat ADCY10 ELISA Kit (orb1197443)

|                            |                                                      |
|----------------------------|------------------------------------------------------|
| <b>Description</b>         | Rat Adenylate Cyclase 10, Soluble (ADCY10) ELISA Kit |
| <b>Reactivity</b>          | Rat                                                  |
| <b>Range</b>               | 0.312-20ng/mL                                        |
| <b>Tested Applications</b> | ELISA                                                |
| <b>Note</b>                | For research use only                                |
| <b>Application notes</b>   | sandwich                                             |
| <b>Sample Types</b>        | tissue homogenates or other biological fluids.       |
| <b>Uniprot ID</b>          | <a href="#">Q9Z286</a>                               |
| <b>Sensitivity</b>         | 0.125ng/mL                                           |
| <b>Entrez</b>              | <a href="#">59320</a>                                |
| <b>Expiration Date</b>     | Please enquire.                                      |