

Product Datasheet

Mouse UQCRC2 ELISA Kit (orb1197477)

Description	Mouse Ubiquinol Cytochrome C Reductase Core Protein II (UQCRC2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q9DB77
Sensitivity	0.063ng/mL
Entrez	67003
Expiration Date	Please enquire.