

Product Datasheet

Mouse UCN2 ELISA Kit (orb1197481)

Description	Mouse Urocortin 2 (UCN2) ELISA Kit
Reactivity	Mouse
Range	1.56-100ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q99ML8
Sensitivity	0.52ng/mL
Entrez	171530
Expiration Date	Please enquire.