

Product Datasheet

Mouse TREM2 ELISA Kit (orb1197497)

Description	Mouse Triggering Receptor Expressed On Myeloid Cells 2 (TREM2) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q99NH8
Sensitivity	23.8pg/mL
Entrez	83433
Expiration Date	Please enquire.