

Product Datasheet

Mouse TPSb2 ELISA Kit (orb1197500)

Description	Mouse Tryptase Beta 2 (TPSb2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P21845
Sensitivity	0.054ng/mL
Expiration Date	Please enquire.