

Product Datasheet

Mouse TNNT2 ELISA Kit (orb1197510)

Description	Mouse Troponin T Type 2, Cardiac (TNNT2) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P50752
Sensitivity	5.3pg/mL
Entrez	21956
Expiration Date	Please enquire.