

## Product Datasheet

### Mouse STC1 ELISA Kit (orb1197534)

|                            |                                           |
|----------------------------|-------------------------------------------|
| <b>Description</b>         | Mouse Stanniocalcin 1 (STC1)<br>ELISA Kit |
| <b>Reactivity</b>          | Mouse                                     |
| <b>Range</b>               | 78.1-5000pg/mL                            |
| <b>Tested Applications</b> | ELISA                                     |
| <b>Note</b>                | For research use only                     |
| <b>Application notes</b>   | sandwich                                  |
| <b>Sample Types</b>        | serum, plasma or other biological fluids. |
| <b>Uniprot ID</b>          | <b>O55183</b>                             |
| <b>Sensitivity</b>         | 29.3pg/mL                                 |
| <b>Entrez</b>              | <b>20855</b>                              |
| <b>Expiration Date</b>     | Please enquire.                           |