

Product Datasheet

Mouse SRSF4 ELISA Kit (orb1197537)

Description	Mouse Serine/Arginine Rich Splicing Factor 4 (SRSF4) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q542V3
Sensitivity	0.043ng/mL
Entrez	57317
Expiration Date	Please enquire.