

## Product Datasheet

### Mouse Slit3 ELISA Kit (orb1197554)

|                            |                                                                                                      |
|----------------------------|------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse Slit Homolog 3 (Slit3) ELISA Kit                                                               |
| <b>Reactivity</b>          | Mouse                                                                                                |
| <b>Range</b>               | 31.2-2000pg/mL                                                                                       |
| <b>Tested Applications</b> | ELISA                                                                                                |
| <b>Note</b>                | For research use only                                                                                |
| <b>Application notes</b>   | sandwich                                                                                             |
| <b>Sample Types</b>        | serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids. |
| <b>Uniprot ID</b>          | <b>Q9WVB4</b>                                                                                        |
| <b>Sensitivity</b>         | 14.5pg/mL                                                                                            |
| <b>Entrez</b>              | <b>20564</b>                                                                                         |
| <b>Expiration Date</b>     | Please enquire.                                                                                      |