

Product Datasheet

Mouse Slit1 ELISA Kit (orb1197555)

Description	Mouse Slit Homolog 1 (Slit1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q80TR4
Sensitivity	0.057ng/mL
Entrez	20562
Expiration Date	Please enquire.