

Product Datasheet

Mouse SFRP5 ELISA Kit (orb1197570)

Description	Mouse Secreted Frizzled Related Protein 5 (SFRP5) ELISA Kit
Reactivity	Mouse
Range	0.25-16ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q9WU66
Sensitivity	0.09ng/mL
Entrez	54612
Expiration Date	Please enquire.