

Product Datasheet

Mouse S100A4 ELISA Kit (orb1197589)

Description	Mouse S100 Calcium Binding Protein A4 (S100A4) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P07091
Sensitivity	14.2pg/mL
Entrez	20198
Expiration Date	Please enquire.