

Product Datasheet

Mouse ROS1 ELISA Kit (orb1197593)

Description	Mouse C-Ros Oncogene 1, Receptor Tyrosine Kinase (ROS1) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q78DX7
Sensitivity	33pg/mL
Entrez	19886
Expiration Date	Please enquire.