

Product Datasheet

Mouse NRG4 ELISA Kit (orb1197665)

Description	Mouse Neuregulin 4 (NRG4) ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q9WTX4
Sensitivity	0.128ng/mL
Entrez	83961
Expiration Date	Please enquire.