

Product Datasheet

Mouse MUSK ELISA Kit (orb1197691)

Description	Mouse Muscle, Skeletal, Receptor Tyrosine Kinase (MUSK) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q61006
Sensitivity	0.061ng/mL
Entrez	18198
Expiration Date	Please enquire.