

Product Datasheet

Mouse LGI3 ELISA Kit (orb1197737)

Description	Mouse Leucine Rich Repeat LGI Family, Member 3 (LGI3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Tested Applications	ELISA
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q8K406
Sensitivity	0.041ng/mL
Entrez	213469
Expiration Date	Please enquire.